## NBF Prep Worksheet (1 L)

Goal: prepare ~4% formaldehyde (10% NBF) from 37% stock.

civi = c2v2.

C1=37, C2=4, V2=1000 mL  $\rightarrow$  V1 =  $(4/37)*1000 \sim 108$  mL formaldehyde stock.

Add phosphate buffer salts to pH ~7.0-7.2; bring to volume with distilled water.

Record: date, operator, lot numbers (formaldehyde, buffer), final pH, label container.

Safety: work in hood; PPE; secondary containment; spill kit.

## Dilution Practice

- 1) How much 37% stock for 500 mL of 10% RBF? V1 ~ 54 mL.
- 2) If final pH reads 6.6: adjust buffer, verify meter, recheck after mixing.
- 3) List three QC items to log for each batch.